

ultra filter ultra fast



with Amicon® Ultra filters

Achieve the fastest and most efficient concentration or buffer exchange for even the most sensitive downstream applications with Amicon® Ultra filters featuring the following benefits:

- **High sample recovery:** 90% recovery from even the most dilute samples thanks to the low analyte binding of the used Ultracel® membranes
- **Reliability:** a true dead stop provides consistent recovery volumes without the risk of overspinning your precious samples to dryness
- **Less clogging, faster filtering:** the dual vertical membrane design is aligned with the filtrate rather than perpendicular to it
- **Quick processing time:** 10–30 minute spin times
- **Maximum sample retrieval:** the integrated reverse spin feature allows for optimal recovery, especially important with dilute and small volume samples

To learn more, please contact us or visit:

SigmaAldrich.com/amicon



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada



Millipore®

Preparation, Separation, Filtration & Monitoring Products